



<b>Department:</b>	Laboratory and Blood Bank (Haematology)		
<b>Document:</b>	Internal Policy and Procedures		
<b>Title:</b>	Giemsa Staining		
<b>Applies To:</b>	All Laboratory Staff		
<b>Preparation Date:</b>	January 06, 2025	<b>Index No:</b>	LB-IPP-044
<b>Approval Date:</b>	January 20, 2025	<b>Version :</b>	2
<b>Effective Date:</b>	February 20, 2025	<b>Replacement No.:</b>	LB-IPP-083 (1)
<b>Review Date:</b>	February 20, 2028	<b>No. of Pages:</b>	02

## 1. PURPOSE:

- 1.1 To establish system and set responsibilities for processing Giemsa stain

## 2. DEFINITIONS:

- 2.1 **Giemsa stain** is used in cytogenetics and for the histopathological diagnosis of malaria and other parasites.

## 3. POLICY:

- 3.1 Place slides into the working Giemsa stain (2.5%) for 45-60 minutes.
- 3.2 The smears could be stained for shorter times in more concentrated stains.

## 4. PROCEDURE:

### 4.1 Preparation of the specimen:

- 4.1.1 Treat specimens containing mucus (sputum mainly) with a mucolytic agent (e.g., Sputalysin) by mixing the agent in a 1:1 ratio with the specimen.
- 4.1.2 Centrifuge at 500 ~ g for 5 minutes.
- 4.1.3 Discard supernatant.
- 4.1.4 Re-suspend the remaining sediment.
- 4.1.5 Place drops of the sediment on the centre of a microscopic slide.
- 4.1.6 Spread the drop on the slide with a pipette to evenly distribute on the slide.
- 4.1.7 Air dry and stain following the required procedure.

### 4.2 Staining:

- 4.2.1 Prepare fresh working Giemsa stain in a staining jar.
- 4.2.2 Place slides into the working Giemsa stain (2.5%) for 45-60 minutes.
- 4.2.3 Remove thin smear slides and rinse by dipping 3-4 times in the Giemsa buffer. Thick smears should be left in buffer for 5 minutes.
- 4.2.4 Dry the slides upright in a rack.
- 4.2.5 As alternates to these 45-60 minutes in 2.5% Giemsa stain, the smears could be stained for shorter times in more concentrated stains; one alternate is 10 minutes in 10% Giemsa; the shorter stains yield faster results, but use more stain and might be of less predictable quality.

## 5. MATERIAL AND EQUIPMENT:

- 5.1 Giemsa stain reagent
- 5.2 Light microscope with an oil immersion lens 100X



## 6. RESPONSIBILITIES:

- 6.1 Assigned Haematology Technician
- 6.2 Clinical Pathology Specialist/ Consultant.







## 7. APPENDICES:

N/A

## 8. REFERENCES:

- 8.1 Procedure Manual, Toronto Medical laboratories / Mount Sinai Hospital department of microbiology.
- 8.2 Bailey & Scott's Diagnostic Microbiology. Feingold & Baron; 12th. Ed.2007, C.V. Mosby Co. p. 301.
- 8.3 Clinical Microbiology Procedures Handbook, American Society of Microbiology, Washington DC, 2005.

## 9. APPROVALS:

	Name	Title	Signature	Date
<b>Prepared by:</b>	Dr. Fatma Hassan Ahmed	Clinical Pathologist		January 06, 2025
<b>Reviewed by:</b>	Dr. Kawther M. Abdou	Consultant & Lab. Medical Director		January 08, 2025
<b>Reviewed by:</b>	Ms. Noora Melfi Alanizi	Laboratory & Blood Bank Director		January 09, 2025
<b>Reviewed by:</b>	Mr. Abdulelah Ayed Al Mutairi	QM&PS Director		January 12, 2025
<b>Reviewed by:</b>	Dr. Tamer Mohamed Naguib	Medical Director		January 12, 2025
<b>Approved by:</b>	Mr. Fahad Hazam Alshammari	Hospital Director		January 20, 2025